

PATIENT LABEL

Bravo® pH Patient Diary

				11	ļ.	=4	Δ	0	0					
Time of Event	Circle am or pm	W. W.	Mear	Swings	Supine Fra	Per Pain	H. H.	Och Com	¢/		Commer	its		
120	am pm		1											_
:	am pm													
:	am pm			ij										
:	am pm													
:	am pm													
:	am pm													
:	am pm													
:	am pm			ij.										
:	am pm													
:	am pm													
:	am pm			Ü										
:	am pm													_
:	am pm													_
:	am pm			j.										_
:	am pm													_
:	am pm		7											_
:	am pm													

Always refer to the time displayed on the recorder to write down the Time of Event above.

Recording event on the Bravo pH recorder

\$78 *	Press any of the three Symptom Buttons (chest pain, regurgitation, hearth urn) once for each individual event.
T 1	Press the Meal Button once at the start of the meal and once at the end of the meal. In between the two presses, the button LED will blink until the button is pressed to mark the end the event.
 4	Press the Supine Button once when lying down for bed and once at the end of the lying period. In between the two presses, the button LED will blink until the button is pressed to mark the end the event.

Notes

All button functions are active only when backlight is ON. If OFF, pressing any of the enabled recorder buttons will first turn backlight ON and no function will be activated. Pressing the desired button when the backlight is ON will activate the desired function.